ELEARNING VITEK 2
MÓDULO 4

ANÁLISIS DE RESULTADOS

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
1 - 20 of 24 items

User Name Laboratory Administrator (LabAdmin)

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Organism: Staph.saprophyticus

AES Findings: Consistent

Phenotypes Selected for Review: OXAZOLIDINONE RESISTANT

Card Comments:

Advanced Reporting Tool Comments:

AST-GP67 !

Antibiotic	MIC	INT	Antibiotic	MIC	INT	Antibiotic	MIC	INT
<input type="checkbox"/> Beta-Lactamase	POS	+	<input type="checkbox"/> Ciprofloxacin	≤0.5	S	<input type="checkbox"/> Vancomycin	1	S
<input type="checkbox"/> Cefoxitin Screen	NEG	-	<input type="checkbox"/> Levofloxacin	0.5	S	<input type="checkbox"/> Tetracycline	≤1	S
<input type="checkbox"/> Benzylpenicillin	0.12	R	<input type="checkbox"/> Moxifloxacin	≤0.25	S	<input type="checkbox"/> Tigecycline	≤0.12	S
<input type="checkbox"/> Ampicillin			<input type="checkbox"/> Inducible Clindamycin Resistance	NEG	-	<input type="checkbox"/> Nitrofurantoin	≤16	S
<input type="checkbox"/> Oxacillin	0.5	S	<input type="checkbox"/> Erythromycin	≤0.25	S	<input type="checkbox"/> Rifampicin	≤0.5	S
<input type="checkbox"/> Gentamicin High Level (synergy)			<input type="checkbox"/> Clindamycin	≤0.25	S	<input type="checkbox"/> Trimethoprim/ Sulfamethoxazole	≤10	S
<input type="checkbox"/> Streptomycin High Level (synergy)			<input type="checkbox"/> Quinupristin/ Dalfopristin	0.5	S			
<input type="checkbox"/> Gentamicin	≤0.5	S	<input checked="" type="checkbox"/> Linezolid	≥8				

Patient Name: Schleicher, Butch

Analysis Status: 14.63 hr - Final

AST Offline Tests: Beta-Lactamase +

Analysis Messages: The following antibiotic(s) are not claimed: Ampicillin, Gentamicin High Level (synergy), Streptomycin High Level (synergy).

McFarland: 0,56

Setup Tech: Laboratory Supervisor (LabSuper)

Organism Quantity: Rto >100.000 ufc

BP Infection Site:

Supplemental Tests:

Contraindicating Tests: Staph.saprophyticus MBdG (14)

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


B1. Enterobacteriaceae

Antimicrobial Agent Organism	Ampicillin	Amoxicillin-clavulanate	Ampicillin-sulbactam	Piperacillin	Ticarcillin	Cephalosporin I: Cefazolin, Cephalothin	Cephameycins: Cefoxitin, Cefotetan	Cephalosporin II: Cefuroxime	Imipenem	Tetracyclines	Tigecycline	Nitrofurantoin	Polymyxin B Colistin	Aminoglycosides
<i>Citrobacter freundii</i>	R	R	R			R	R	R						
<i>Citrobacter koseri</i>	R			R	R									
<i>Enterobacter aerogenes</i>	R	R	R			R	R	R						
<i>Enterobacter cloacae</i> complex	R	R	R			R	R	R						
<i>Escherichia coli</i>	There is no intrinsic resistance to β -lactams in this organism.													
<i>Escherichia hermannii</i>	R				R									
<i>Hafnia alvei</i>	R	R	R			R	R							
<i>Klebsiella pneumoniae</i>	R				R									
<i>Morganella morganii</i>	R	R				R		R	*		R	R	R	
<i>Proteus mirabilis</i>	There is no intrinsic resistance to penicillins and cephalosporins in this organism.									*	R	R	R	
<i>Proteus penneri</i>	R					R		R	*	R	R	R	R	
<i>Proteus vulgaris</i>	R					R		R	*	R	R	R	R	
<i>Providencia rettgeri</i>	R	R				R			*	R	R	R	R	
<i>Providencia stuartii</i>	R	R				R			*	R	R	R	R	†
<i>Salmonella</i> and <i>Shigella</i> spp.	There is no intrinsic resistance to β -lactams in these organisms; refer to WARNING below for reporting.													
<i>Serratia marcescens</i>	R	R	R			R	R	R				R	R	
<i>Yersinia enterocolitica</i>	R	R			R	R								

WARNING: For *Salmonella* spp. and *Shigella* spp., aminoglycosides, first- and second-generation cephalosporins, and cephamycins may appear active *in vitro*, but are not effective clinically and should not be reported as susceptible.

* *Proteus* species, *Providencia* species, and *Morganella* species may have elevated minimal inhibitory concentrations to imipenem by mechanisms other than by production of carbapenemases. Isolates that test as susceptible should be reported as susceptible.

† *Providencia stuartii* should be considered resistant to gentamicin, netilmicin, and tobramycin but not intrinsically resistant to amikacin.

NOTE 1: Cephalosporins III, cefepime, aztreonam, ticarcillin-clavulanate, piperacillin-tazobactam, and the carbapenems are not listed, because there is no intrinsic resistance in *Enterobacteriaceae*.

NOTE 2: *Enterobacteriaceae* are also intrinsically resistant to clindamycin, daptomycin, fusidic acid, glycopeptides (vancomycin, telavancin), lipoglycopeptides (oritavancin, telavancin), linezolid, tedizolid, quinupristin-dalfopristin, rifampin, and macrolides (erythromycin, clarithromycin, and azithromycin). However, there are some exceptions with macrolides (ie, *Salmonella* and *Shigella* spp. with azithromycin).

Laboratory Standards Institute antimicrobial susceptibility testing standards M02-A12, M07-A10, and M11-A8.

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* *Acinetobacter baumannii/calcoaceticus* may appear to be susceptible to ampicillin-sulbactam due to the activity of sulbactam with this species.

† *Stenotrophomonas maltophilia* is intrinsically resistant to tetracycline but not to doxycycline, minocycline, or tigecycline.

NOTE: These nonfermentative gram-negative bacteria are also intrinsically resistant to penicillin (ie, benzylpenicillin), cephalosporin I (cephalothin, cefazolin), cephalosporin II (cefuroxime), cephamycins (cefoxitin, cefotetan), clindamycin, daptomycin, fusidic acid, glycopeptides (vancomycin, teicoplanin), linezolid, macrolides (erythromycin, azithromycin, clarithromycin), quinupristin-dalfopristin, and rifampin.

Laboratory Standards Institute antimicrobial susceptibility testing standards M02-A12, M07-A10, and M11-A8.

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NOTE 1: These gram-positive bacteria are also intrinsically resistant to aztreonam, polymyxin B/colistin, and nalidixic acid.

NOTE 2: Oxacillin-resistant *S. aureus* and coagulase-negative staphylococci (methicillin-resistant staphylococci [MRS]) are considered resistant to other β -lactam agents, ie, penicillins, β -lactam/ β -lactamase inhibitor combinations, cepheims (with the exception of the cephalosporins with anti-MRSA [methicillin-resistant *S. aureus*] activity), and carbapenems. This is because most cases of documented MRS infections have responded poorly to β -lactam therapy, or because convincing clinical data that document clinical efficacy for those agents have not been presented.

Laboratory Standards Institute antimicrobial susceptibility testing standards M02-A12, M07-A10, and M11-A8.

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* **Warning:** For *Enterococcus* spp., cephalosporins, aminoglycosides (except for high-level resistance testing), clindamycin, and trimethoprim-sulfamethoxazole may appear active *in vitro*, but are not effective clinically and should not be reported as susceptible.

NOTE: These gram-positive bacteria are also intrinsically resistant to aztreonam, polymyxin B/colistin, and nalidixic acid.

Laboratory Standards Institute antimicrobial susceptibility testing standards M02-A12, M07-A10, and M11-A8.

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